



<b>DE</b>	Diagnosewaage Gebrauchsanweisung .....	3
<b>EN</b>	Diagnostic bathroom scale Instructions for use .....	17
<b>FR</b>	Pèse-personne impédancemètre Mode d'emploi .....	30
<b>ES</b>	Báscula diagnóstica Instrucciones de uso.....	43
<b>IT</b>	Bilancia diagnostica Istruzioni per l'uso.....	57
<b>TR</b>	Diyanoz terazisi Kullanım kılavuzu .....	70
<b>RU</b>	Диагностические весы Инструкция по применению .....	82
<b>PL</b>	Waga diagnostyczna Instrukcja obsługi .....	96



**A**

**B**

**C**

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**Unfold page 3 before reading the instructions for use.**

**Read these instructions for use carefully. Observe the warnings and safety notes. Keep these instructions for use for future reference. Make the instructions for use accessible to other users. If the device is passed on, provide the instructions for use to the next user as well.**

## **INCLUDED IN DELIVERY**

- 1 x BF 950 diagnostic bathroom scale
- 1 x 4 x 1.5 V AAA batteries
- 1 x Feet for carpet
- 1 x These instructions for use
- 1 x Warranty leaflet
- 2 x Quick guide

## **Contents**

1. IMPORTANT SAFETY NOTES .....	18
2. SIGNS AND SYMBOLS .....	18
3. DEVICE DESCRIPTION .....	19
4. INFORMATION .....	19
5. INITIAL USE OF THE SCALE WITH THE APP .....	20
6. INITIAL USE OF THE SCALE WITHOUT THE APP .....	21
7. TAKING MEASUREMENTS .....	22
8. EVALUATING THE RESULTS .....	23
9. ADDITIONAL INFORMATION .....	26
10. CLEANING AND MAINTENANCE .....	26
11. DISPOSAL .....	27
12. WHAT IF THERE ARE PROBLEMS? .....	27
13. TECHNICAL SPECIFICATIONS .....	28
14. WARRANTY .....	29

# 1. IMPORTANT SAFETY NOTES

## WARNING

- The device is only intended for weighing humans and for recording your personal fitness data. The device is only intended for private use and not for medical or commercial purposes. The device must only be used for the purpose for which it was designed and in the manner specified in the instructions for use.
- Danger of suffocation! Keep packaging materials away from children.
- **The scale may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.** 
- When pregnant, only use the scale in pregnancy mode (see the “Pregnancy mode” chapter).
- Danger of tipping! Never place the scale on an uneven surface. Never step onto the outer edge of the scale on just one side (see Figure ).
- Danger of slipping! Do not step onto the scale with wet feet (see Figure ). Do not step onto the scale when the surface is wet. Do not step onto the scale wearing socks.

## Safety notes on handling batteries

- If your skin or eyes come into contact with battery fluid, rinse the affected area with water and seek medical assistance.
-  **Choking hazard!** Small children may swallow and choke on batteries. Store batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
-  **Risk of explosion!** Do not throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, open or crush the batteries..

## 2. SIGNS AND SYMBOLS

The following symbols are used on the device, in these instructions for use, on the packaging and on the type plate for the device:

	<b>Warning</b> Warning notice indicating a risk of injury or damage to health.
	<b>Important</b> Safety note indicating possible damage to the device/accessory.
	<b>Note</b> Note on important information
	<b>CE labelling</b> This product satisfies the requirements of the applicable European and national directives.
	The products demonstrably meet the requirements of the Technical Regulations of the EAEU.

	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE		
	United Kingdom Conformity Assessed Mark		
	Manufacturer		Dispose of packaging in an environmentally friendly manner
	Read the instructions		<b>Danger of slipping:</b> do not step on the scale with wet feet.
	Place the scale on an even surface. No carpet.		<b>Danger of tipping:</b> position yourself in the centre of the weighing surface.

### 3. DEVICE DESCRIPTION

- 1** Down button
- 2** SET button
- 3** Up button
- 4** Electrodes (ITO)
- 5** LED motivator:
- 6** Display
- 7** Unit/Delete/Reset button
- 8** Battery compartment
- 9** *Bluetooth*<sup>®</sup> symbol, height in cm/ft, weight in kg/lb/st
- 10** **BMI** = Interpretation of weight  
**BF** = Body fat in %  
 = Body water in %  
 = Muscle percentage in %  
 = Bone mineral percentage in %  
**BMR** = Basal metabolic rate in kcal  
**AMR** = Active metabolic rate in kcal  
**°C/°F** = Temperature in °C/°F
- 11** Initials of the user
- 12**  = Tare function  
 = Pregnancy mode  
 = Female  
 = Male
- 13** Date (Y/M/D)
- 14** Time (12/24 h)

### 4. INFORMATION

#### The measurement principle

This scale works on the principle of B.I.A. (bioelectrical impedance analysis). This involves the calculation of body content in a matter seconds, using a current that cannot be felt, is completely harmless and does not pose any risk.

When this measurement of the electrical resistance (impedance) is considered alongside constants and/or individual values (age, height, gender, activity level), the body fat percentage and other variables in the body

can be calculated. Muscle tissue and water conduct electricity well and therefore have a lower resistance. In contrast, bones and fat tissue have low conductivity as the fat cells and bones are barely able to conduct the current due to very high levels of resistance.

Please note that the values calculated by the diagnostic bathroom scale only represent an approximation of the actual medical analytical values for the body. Only a specialist doctor can precisely establish the body fat, body water, muscle percentage and bone structure using medical methods (e.g. computer tomography).

### General tips

- As far as possible, weigh yourself at the same time everyday (ideally in the morning), after having been to the toilet, on an empty stomach and without clothing, in order to obtain results which can be compared.
- Important point regarding the measurement: the calculation of body fat may only be made when barefoot; the soles of the feet may be lightly dampened where appropriate.
- There must be no skin contact between your feet, legs, calves or thighs. Otherwise, the measurement cannot be taken correctly.
- Stand up straight and still during measurement.
- Wait a few hours after any physical exertion to which your body is unaccustomed.
- Wait approx. 15 minutes after getting up so that the water stored in your body can disperse.
- Remember that only the long-term trend is important. Short-term changes in weight within a few days are normally caused by loss of fluids; however, body water plays an important role in terms of well-being.
- The total of all analytical values comes to more than 100% because body water must not be added; body water is already included in the other analytical values (e.g. body fat, muscle percentage).

### Limitations

When calculating body fat and the other values, unusual and implausible results may occur in the case of:

- Children under approx. 5 years of age
- Performance athletes and bodybuilders
- Persons with a fever, persons undergoing dialysis, persons with oedema symptoms and persons suffering from osteoporosis
- Persons taking cardiovascular medication (affecting the heart and vascular system)
- Persons taking vasodilative or vasoconstrictive medication
- Persons with significant anatomical abnormalities in their legs compared with their overall body height (leg length significantly shorter or longer)

## 5. INITIAL USE OF THE SCALE WITH THE APP

To calculate your body fat percentage and additional body values, you must save the personal user data on the scale.

The scale has 8 user memory spaces which allows, for example, you and your family members to save personal settings.

It is also possible to activate users via other mobile end devices on which the “beurer HealthManager Pro“ / “beurer HealthManager” app has been installed or by changing the user profile in the app.

Keep the smartphone close to the scale in order to maintain an active *Bluetooth*<sup>®</sup> connection during initial use.

- 1. Activate *Bluetooth*<sup>®</sup> in your smartphone’s settings.**
- 2. Insert the batteries and place the scale on a solid surface (“USE APP” will flash on the scale’s LCD).**
- 3. Download the free “beurer HealthManager Pro“ / “beurer HealthManager” app from the Apple App Store (iOS) or from Google Play (Android).**
- 4. Start the app and follow the instructions.**

## 5. Enter the BF950-related data.

The following settings must be set or entered in the “beurer HealthManager Pro“ / “beurer HealthManager” app:

User data	Set values												
Gender	Male (♂); female (♀) If user is female: pregnant (🤰) yes/no												
Initials	Max. 3-digit (A-Z, 0-9)												
Height	40 to 250 cm (1' 3.7" to 8' 2.4")												
Date of birth	Year, month, day If user <2 years old: Tare function (👶) yes/no												
Activity level	The medium-term and long-term perspective is decisive when selecting the activity level: <table border="1"><thead><tr><th colspan="2">Physical activity</th></tr></thead><tbody><tr><td>1</td><td>None.</td></tr><tr><td>2</td><td>Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).</td></tr><tr><td>3</td><td>Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.</td></tr><tr><td>4</td><td>High: physical exertion, at least 4 to 6 times a week for 30 minutes each.</td></tr><tr><td>5</td><td>Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.</td></tr></tbody></table>	Physical activity		1	None.	2	Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).	3	Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.	4	High: physical exertion, at least 4 to 6 times a week for 30 minutes each.	5	Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.
Physical activity													
1	None.												
2	Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).												
3	Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.												
4	High: physical exertion, at least 4 to 6 times a week for 30 minutes each.												
5	Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.												
Target weight	5–180 kg or no target. If “pregnant” has been selected under gender, no target weight can be specified.												

The app assigns the next free scale user memory (P01 – P08) to you.

## 6. Carrying out the initial measurement

For automatic personal recognition, the first measurement must be assigned to your personal user data. The app will prompt you to do this. Step onto the scale, ideally with bare feet, and ensure that you are standing still with equal weight distribution and with both feet on the electrodes. If you have selected the tare function, you must carry out weighing twice, once with and once without the child.

## 7. Now continue with chapter “Taking measurements”.

# 6. INITIAL USE OF THE SCALE WITHOUT THE APP

We recommend using the app for the initial use and implementing all settings. See chapter 6. “Initial use of the scale with the app”. Alternatively, you can configure settings (time, date, creating a user) on the scale itself.

## Setting the time and date

1. Remove the batteries from the protective packaging and insert the batteries into the scale, ensuring that the battery polarity is correct. “USE APP” flashes on the scale LCD.
2. Place the scale on an even, hard surface; a hard surface is crucial for accurate measurements.
3. Press the “SET” button.
4. Use the “∨ ^” buttons to set the date (year/month/day) and time (hour/minute) and confirm with the “SET” button.
5. To change the weight unit (kg/lb/st), briefly press the Unit/Delete/Reset button when the scale is switched on.

## Creating a user

1. Briefly apply pressure to the scale twice until “0.0” is displayed.
2. Press the “SET” button. “P-01” is displayed.
3. Use the “∨ ^” buttons to select the user memory space you wish to use (P-01 to P-08). Confirm using the “SET” button.

- Use the “√ ^” to configure your user data and confirm with the “SET” button. You can find a detailed description of all the user data to be configured in the table in chapter “Initial use of the scale with the app”.
- A PIN will be displayed at the end of your settings. This will be required for logging into the app later (see “Logging into the app at a later time” below); otherwise, it is not needed. “0.0” with initials then appears.
- For automatic personal recognition, the first measurement must be assigned to your personal user data. Step onto the scale, ideally with bare feet, and ensure that you are standing still with equal weight distribution and with both feet on the electrodes.

### Logging into the app at a later time

If you have carried out the initial use of the scale without using the app, you can transfer or create your user on the app as follows.

- Download the free “beurer HealthManager Pro“ / “beurer HealthManager” app from the Apple App Store (iOS) or from Google Play (Android).
- Start the app and follow the instructions.
- Select the existing scale user within the BF950 settings in the app and enter the PIN shown on the scale display. Briefly apply pressure to the scale twice until “0.0” is displayed. Press the “SET” button and select your user with the “SET” button.

## 7. TAKING MEASUREMENTS

### Diagnostic functions


Step onto the scale with bare feet and ensure that you are standing still with equal weight distribution and with both legs on the electrodes.

**Note:** The measurement result will be incorrect if the measurement is taken with socks on.

The time and temperature are shown first for 3 seconds, then the weight.


If a user has been assigned, their initials are displayed. During the BIA analysis, “0000” is displayed. All of the following values are displayed and stored. All values are displayed and saved consecutively.

**Note:** If the user is not recognised, then the user assignment (initials) is missing and only the weight can be displayed. The measured values will not be saved. This happens if your measured weight differs from the previous measurement by +/- 3 kg or if another user has similar values (+/- 3 kg) to yours.


In this case, repeat the measurement by turning on the scale and pressing the “^” button to select your user memory (P01 – P08). Repeat the measurement.


### Weight

Step onto the scale wearing shoes and ensure that you are standing still with equal weight distribution and with both legs on the scale.

The time and temperature are shown first for 3 seconds, then the weight.

If a user was assigned, your initials will be displayed and the weight will be saved. "--" will briefly be displayed on the LCD. A BIA analysis is not possible.

## Tare function

The tare function provides an easy method of weighing small children (under 2 years old) using tare weighing. To do this, you must create a user with user data for the infant in the “beurer HealthManager Pro“ / “beurer HealthManager” app. Alternatively, you can create a user with user data on the scale itself (see “Creating a user”).

In the tare function, automatic user assignment for this user is switched off. To use the tare function, proceed as follows for each measurement:

1. Briefly apply pressure to the scale.
2. Use the “^” button to select the user memory for the infant and wait 3 seconds. The tare function symbol  appears on the display.
3. Now step on the scale. You can choose whether to hold the infant in your arms on the first or second measurement.
4. After the first measurement, the first weight is briefly displayed. Step off the scale; “0.0” appears and the scale is ready for the second measurement. Both values appear on the display first; the display then switches to the difference in weight, which is also saved.

## Pregnancy mode

Pregnancy mode allows pregnant women to monitor and compare their weight in the app with greater precision during this special time.

To do so, the pregnant woman must create a user in the “beurer HealthManager Pro“ / “beurer HealthManager” app and select “pregnant” in the user data. If you select “pregnant” in the user data, you cannot select a target weight.

Alternatively, you can create a user with user data on the scale itself (see “Creating a user”).

When pregnancy mode is activated, no current flows through the body during weighing and no analytical values can be determined. Only the weight with initials is displayed and the pregnancy symbol  flashes.

To use pregnancy mode, proceed as follows:

1. Briefly apply pressure to the scale.
2. Use the “^” button to select the user memory for the pregnant woman. The pregnancy symbol  appears on the display.
3. Step onto the scale. The weight appears on the display and the pregnancy symbol  flashes briefly.

## 8. EVALUATING THE RESULTS

### Target weight LED motivator

The “target weight” function must be activated in the app settings (or using the buttons on the scale) and a corresponding target weight must have been stored in the respective user data. If “pregnant” has been specified in the user data, no target weight can be stored.

The LED motivator provides coloured LED feedback after weighing to indicate whether the target weight has been achieved.

There are various scenarios:

Weight maintained:	The target weight is less than +/- 1 kg different to the start value. <ul style="list-style-type: none"><li>• All 9 LEDs light up in green if the next measurements are within +/-1 kg of the start value.</li><li>• The first LED lights up in red if the measurements are outside this range.</li></ul>
--------------------	---

Weight changed (lost or gained):	<p>The target weight is more than +/- 1 kg different to the start value.</p> <ul style="list-style-type: none"> <li>• If the current measurement is closer to the target weight, this is indicated in green. The more LEDs light up green, the closer you are to your target weight (linear distribution). When there are 9 green LEDs, you have achieved your target.</li> <li>• If the current measurement is within +/- 0.2 kg respectively of the start value (moving away from the target value), the first LED lights up in yellow.</li> <li>• If the current measurement is more than 0.2 kg from the start value (moving away from the target value), the first LED lights up in red.</li> </ul>
----------------------------------	--

### Body mass index

The body mass index (BMI) is a figure that is commonly used to evaluate body weight. The figure is calculated from body weight and height. The formula is:  $\text{body mass index} = \text{body weight} \div \text{height}^2$ . The measurement unit for BMI is [kg/m<sup>2</sup>]. According to the BMI, weight is classified for adults (20 years and over) using the following values:

Category		BMI
Underweight	Severely underweight	<16
	Moderately underweight	16-16.9
	Slightly underweight	17-18.4
Normal weight		18.5-24.9
Overweight	Pre-obese	25-29.9
Obese (overweight)	Class I obesity	30-34.9
	Class II obesity	35-39.9
	Class III obesity	≥40

### Body fat percentage

Values will often be lower for fitness enthusiasts. Depending on the sporting activities performed, the training intensity and the physical constitution, the results may still be below the stated standard values. Please note, however, that there may be a risk to health in the case of extremely low values.

The subsequent body fat percentages serve as a guideline (for further information ask your doctor).

#### Male

Age	Low	Normal	High	Very high
10-14	<11 %	11-16 %	16,1-21 %	<21 %
15-19	<12 %	12-17 %	17,1-22 %	<22 %
20-29	<13 %	13-18 %	18,1-23 %	<23 %
30-39	<14 %	14-19 %	19,1-24 %	<24 %
40-49	<15 %	15-20 %	20,1-25 %	<25 %
50-59	<16 %	16-21 %	21,1-26 %	<26 %
60-69	<17 %	17-22 %	22,1-27 %	<27 %
70-100	<18 %	18-23 %	23,1-28 %	<28 %

#### Female

Age	Low	Normal	High	Very high
10-14	<16 %	16-21 %	21,1-26 %	<26 %
15-19	<17 %	17-22 %	22,1-27 %	<27 %
20-29	<18 %	18-23 %	23,1-28 %	<28 %
30-39	<19 %	19-24 %	24,1-29 %	<29 %
40-49	<20 %	20-25 %	25,1-30 %	<30 %
50-59	<21 %	21-26 %	26,1-31 %	<31 %
60-69	<22 %	22-27 %	27,1-32 %	<32 %
70-100	<23 %	23-28 %	28,1-33 %	<33 %

## Body water content

The body water content in % is normally within the following ranges:

### Male

Age	Poor	Good	Very good
10-100	<50 %	50-65 %	<65 %

### Female

Age	Poor	Good	Very good
10-100	<45 %	45-60 %	<60 %

Body fat contains relatively little water, therefore body water content may be below the standard values in persons with a high body fat percentage. On the other hand, endurance athletes may exceed the standard values due to having low levels of fat and a high muscle percentage.

The body water calculation performed using this scale is not suitable for drawing medical conclusions, for example in the case of age-related water retention. Consult your doctor where necessary. Body water must not be added to the other values, as this is already included in the other analytical values (e.g. body fat, muscle percentage).

## Muscle percentage

The muscle percentage measurement measures the skeletal muscle mass. No organ muscles are taken into account. The muscle percentage is normally within the following ranges:

### Male

Age	Low	Normal	High
10-14	<44 %	44-47 %	>57 %
15-19	<43 %	43-56 %	>56 %
20-29	<42 %	42-54 %	>54 %
30-39	<41 %	41-52 %	>52 %
40-49	<40 %	40-50 %	>50 %
50-59	<39 %	39-48 %	>48 %
60-69	<38 %	38-47 %	>47 %
70-100	<37 %	37-46 %	>46 %

### Female

Age	Low	Normal	High
10-14	<36 %	36-43 %	>43 %
15-19	<35 %	35-41 %	>41 %
20-29	<34 %	34-39 %	>39 %
30-39	<33 %	33-38 %	>38 %
40-49	<31 %	31-36 %	>36 %
50-59	<29 %	29-34 %	>34 %
60-69	<28 %	28-33 %	>33 %
70-100	<27 %	27-32 %	>32 %

## Bone mass

Like the rest of our body, our bones are subject to natural growth, deterioration and ageing processes. Bone mass increases rapidly during childhood and reaches its maximum value between the ages of 30 and 40. As we age, our bone mass then gradually reduces. A healthy diet (particularly one rich in calcium and vitamin D) and regular exercise can slow down this process. The stability of the skeletal structure can also be boosted by targeted muscle build-up.

Please note that this scale shows the mineral content of bone and not any organic substances such as proteins or cells. It is virtually impossible to influence bone mass; yet there may be distinct fluctuations due to certain influential factors (such as weight, height, age, gender). There are no recognised guidelines or recommendations.

## BMR

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest in order to maintain its basic functions (e.g. when lying in bed 24 hours a day). This value is primarily dependent on weight, height and age.

This is the amount of energy that is required by your body under all circumstances and must be re-supplied to the body in the form of food. If your energy intake is below this level in the longer term, this can result in damage to health.

## AMR

The active metabolic rate (AMR) is the amount of energy that the active body consumes each day. An individual's energy requirement increases when the level of physical activity increases; the diagnostic bathroom scale calculates this by means of the activity level (1-5) that has been entered.

If less energy is taken in than is used over a long period of time, the body primarily takes the difference from the fat stores with resulting weight loss. However, if the amount of energy taken in exceeds the calculated active metabolic rate (AMR) for a longer period, the body cannot burn off the excess energy. The excess is stored in the body as fat, leading to weight gain.

### Temporal context of the results

**Note:** Please note that only the long-term trend is important. Brief deviations in weight within a few days are normally caused by loss of fluids.

## 9. ADDITIONAL INFORMATION

### User assignment

It is possible to assign up to 8 scale users for a weight-only measurement (with shoes) and for a diagnosis measurement (bare feet). User assignment is not possible with the tare function.

For a new measurement, the scale assigns the measurement to the user whose most recently saved measurement falls within +/- 3kg of this.

### Storing the assigned measurements on the scale and transferring them to the app.

If the app is not open, the newly assigned measurements are saved on the scale. A total of 30 measurements per user can be saved on the scale. The saved measurements are transferred automatically to the app when you open the app within the Bluetooth range.

### Deleting scale data

- **Delete all**

Briefly apply pressure to the scale; the time and temperature are displayed. Hold down the Unit/Delete/Reset button on the rear of the scale for approx. 6 seconds. "CLr" will appear on the display for several seconds.

**Note:** All saved values and settings are deleted.

- **Delete individual users**

Briefly apply pressure to the scale twice until "0.0" appears. Press the "√ ^" buttons to select your user. When "0.0" and your initials are displayed, hold down the Unit/Delete/Reset button on the rear of the scale for approx. 3 seconds. "dEL" appears on the display.

**Note:** This user is deleted.

### Changing the batteries

Your scale is equipped with a low battery indicator. If you operate the scale with flat batteries, "Lo" will appear on the display and the scale will automatically switch off. In this case, the batteries must be replaced (4 x 1.5 V AAA). All saved measurements and users remain saved.

Connect the scale to the app again. The system time is updated.

## 10. CLEANING AND MAINTENANCE

The device should be cleaned from time to time. This should be done using a damp cloth and, if necessary, a small amount of detergent.

### IMPORTANT

- Never use harsh solvents or cleaning agents!
- Never submerge the device in water!
- Do not clean the device in the dishwasher.

The accuracy of the measurements and service life of the device depend on its careful handling:

- Do not place any objects on the scale when not in use.
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).
- Do not press the buttons violently or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields.

## 11. DISPOSAL

Empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead

Cd = Battery contains cadmium

Hg = Battery contains mercury

For environmental reasons, do not dispose of the device in the household waste at the end of its service life. Dispose of the device at a suitable local collection or recycling point in your country. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Dispose of packaging in an environmentally friendly manner

## 12. WHAT IF THERE ARE PROBLEMS?

If the scale encounters an error during measurement, the following is displayed:

Display	Cause	Solution
<b>Only weight display and initials are missing</b>	Unknown measurement as it is outside the user assignment limit or a unique assignment is not possible.	Briefly apply pressure to the scale and select the user on the scale using the “^” button or repeat user assignment in the app. See chapter 10 “Additional information”.
Err	The proportion of fat is outside the measurable range (less than 4% or greater than 65%).	Please repeat the measurement barefoot, or with the soles of the feet lightly dampened if necessary. If Err is still displayed, you are outside the measurable range.
Old	Maximum weight capacity of 180 kg exceeded.	Loads must be limited to 180 kg.
<b>Incorrect weight is displayed.</b>	No flat, stable surface.	Place the scale on a flat, stable surface. Do not use on carpet. If used on carpet, use the supplied feet for carpet.
<b>Incorrect weight is displayed.</b>	Not standing still.	Stand as still as possible.
<b>Incorrect weight is displayed.</b>	Incorrect zero point of scale.	Wait until the scale switches itself off. Briefly apply pressure to the scale until “0.0” appears. Then repeat the measurement.

Display	Cause	Solution
<b>No Bluetooth® connection</b> (📶 symbol is missing).	Device outside the range.	The minimum range in the open air is approx. 10 m. Walls and ceilings reduce the range. Other radio waves may interfere with the transmission. Therefore, do not position the scale near devices such as micro-waves, induction hobs etc. and, ideally, switch off any other <i>Bluetooth</i> ® devices that are within range.
	Connection set-up takes too long.	Faster <i>Bluetooth</i> ® interval when the scale is switched on.
	No connection to the app.	<ul style="list-style-type: none"> <li>• Close all apps completely (including in the background).</li> <li>• Switch <i>Bluetooth</i>® off and back on again.</li> <li>• Switch the smartphone off and back on again.</li> <li>• In the <i>Bluetooth</i>® settings of the operating system, the scale must not be connected.</li> <li>• Remove the batteries from the scale for a short period and reinsert them.</li> <li>• See FAQs at <a href="http://www.beurer.com">www.beurer.com</a>.</li> </ul>
FULL	User memory space is full. New measurements overwrite the older measurements.	Open the app. The data is transferred automatically.
Lo	The batteries in the scale are empty.	Replace the batteries in the scale.

### 13. TECHNICAL SPECIFICATIONS

Scale dimensions	320 x 320 x 24 mm
Scale batteries	4 x 1.5V type AAA
Measurement range	5 kg to 180 kg
Scale interval d:	The display can be read in increments of 0.1 kg.
Repetition accuracy:	The measuring tolerance for repeated measurements is +/- 0.4 kg (several measurements in direct succession on the same scale with position of scale and person the same as far as possible).
Absolute precision:	In comparison with a calibrated weight, the measured value is +/-1% +0.1 kg. E.g. at 40 kg this corresponds to +/- 0.5 kg; at 100 kg this corresponds to +/-1.1 kg.
Data transfer via <i>Bluetooth</i> ® wireless technology	This product uses <i>Bluetooth</i> ® low energy technology, frequency band 2.400 – 2.480 MHz, transmission power max. -0,82dBm, compatible with <i>Bluetooth</i> ® 4.0 smartphones/tablets
A list of compatible smartphones, plus information and software for the devices can be found at the following link: <a href="http://connect.beurer.com">connect.beurer.com</a>	

System requirements for the “beurer HealthManager Pro“ / “beurer Health- Manager” App	iOS ≥ 12.0, Android™ ≥ 8.0 <i>Bluetooth</i> ® ≥ 4.0
---	--

Subject to technical changes.

We hereby confirm that this product complies with the European RED Directive 2014/53/EU. The CE Declaration of Conformity for this product can be found under:

[www.beurer.com/web/we-landingpages/de/cedeclarationofconformity.php](http://www.beurer.com/web/we-landingpages/de/cedeclarationofconformity.php).

## 14. WARRANTY

Further information on the warranty and warranty conditions can be found in the warranty leaflet supplied.

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The *Bluetooth*<sup>®</sup> word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

Beurer GmbH • Söflinger Str. 218 • 89077 Ulm (Germany)  
[www.beurer.com](http://www.beurer.com) • [www.beurer.healthguide.com](http://www.beurer.healthguide.com)  
[www.beurer-gesundheitsratgeber.com](http://www.beurer-gesundheitsratgeber.com)