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IPL 6500 Long-lasting hair removal

Compact device for long-lasting hair removal - long-lasting results over several months

Compact design - also ideal for treating small body surfaces

Includes cartridge with 100,000 light pulses

Around just 26 min. to treat the entire body¹⁾

3.1 cm² light surface for especially fast treatment

Free beurer MyIPL app guides through the treatment and creates an individual treatment plan - for a professional IPL treatment at home

Clinically tested technology for maximum safety at home

Maximum safety with 2-in-1 skin-type and skin-contact sensor

Up to 50% hair reduction after 3-4 treatments

State-of-the-art IPL light technology used by dermatologists for professional hair removal. Deactivates the hair roots deep beneath the skin and in doing so decreases the chance of the hair growing back.

With "flash-glide" function for quick and convenient use

Integrated UV filter to protect the skin

5-year warranty


Cannot be used on dark skin

May be ineffective on white or light-coloured hair

Product measurements (LxWxH): 72 x 132 x 44 mm

Product weight:

Device: approx. 224 g / Mains part: approx. 140 g

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 57512 0

Item no.: 575.12

* Dermatological study, 20 women, Institut proDERM Hamburg, study no. 14.0309-11
1) at minimum power level for arms, legs, face, armpits and the bikini line

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