beauty

beurer

**CLINICALLY TESTED TECHNOLOGY**

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works with  
  
 MyIPL  
 Download on the App Store  
 GET IT ON Google play

1 Impuls pro sek.  
 3,1 cm<sup>2</sup>  
 150.000 LICHTIMPULSE  
 LIGHT PULSES  
 3,1 cm<sup>2</sup>  
 1 impulse per sec.

**SIZE compact**  
**long-lasting hair removal**  
**safe skin sensor**

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
**>95%** vertragen die Anwendung sehr gut  
 found their skin tolerated the application very well

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
## IPL 7500 SatinSkin Pro Long-lasting hair removal

- Compact device for long-lasting hair removal - long-lasting results over several months
- Compact design - also ideal for treating small body surfaces
- Includes cartridge with 150,000 light pulses
- Fastest compact device - Around just 20 min. to treat the entire body<sup>1)</sup>
- 3.1 cm<sup>2</sup> light surface for especially fast treatment
- Free beurer MyIPL app guides through the treatment and creates an individual treatment plan - for a professional IPL treatment at home
- Clinically tested technology for maximum safety at home
- Maximum safety with 2-in-1 skin-type and skin-contact sensor
- Up to 50% hair reduction after 3-4 treatments
- State-of-the-art IPL light technology used by dermatologists for professional hair removal. Deactivates the hair roots deep beneath the skin and in doing so decreases the chance of the hair growing back.
- With "flash-glide" function for quick and convenient use
- Integrated UV filter to protect the skin
- 5-year warranty
- Cannot be used on dark skin
- Product measurements (LxWxH): 72 x 132 x 44 mm
- Product weight:  
Device: approx. 224 g / Mains part: approx. 176 g
- May be ineffective on white or light-coloured hair

VE: 2 / Transportkarton: -  
 EAN-Nr.: 4211125 57616 5  
 Art.-Nr.: 576.16

\* Dermatological study, 20 women, Institut proDERM Hamburg, study no. 14.0309-11  
 1) at minimum power level for arms, legs, face, armpits and the bikini line

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