


MG 89  
CompactPower

beurer



## MG 89 MASSAGE GUN CompactPower

### Deep tissue massager



Activation and regeneration through targeted muscle massage



Compact design – perfect for home use, travel or training



3 high-quality massage attachments for personalised usage



Durable thanks to brushless motor

Powerful and soothing trigger point massage

Quiet operation: 46 dB at the first operating level

5 intensity levels:

Intensity level 1: approx. 1800 rpm or 30 Hz

Intensity level 2: approx. 2100 rpm or 35 Hz

Intensity level 3: approx. 2400 rpm or 40 Hz

Intensity level 4: approx. 2700 rpm or 45 Hz

Intensity level 5: approx. 3000 rpm or 50 Hz

Lightweight and handy design for ease of use

Powerful lithium-ion battery with up to 4 hours of battery life

LED lights to display battery and intensity levels

Product measurements (L x W x H): approx. 12,7 x 5,1 x 15,3 cm (without attachments)

Product weight: approx. 431 g (without attachments, without USB-C cable)

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 10434 4

Item no.: 104.34



104.34\_MG89\_2024-10-23\_02\_DS\_EN\_BEU\_Subject to errors and changes