

HK 58
cosy

beurer
wellbeing



HK 58 Cosy Back and neck heat pad




cosy cover

Soft and cosy heat pad with cover



machine washable

Machine washable up to 30 °C



6 temperature settings

6 illuminated temperature settings



Formfitted form for back and neck

Individually adjustable with clip and hook-and-loop fastener

Innovative button control with LED lights

Breathable, cosy and skin-friendly

6 illuminated temperature settings

Automatic switch off after approx. 90 minutes

Rapid heat up

BSS (Beurer Safety System)

Electronic temperature control

Detachable control unit

Machine washable up to 30 °C (heat pad + cover)

OEKO-TEX®

The textiles used for this device meet the stringent human ecological requirements of OEKO-TEX® Standard 100, as verified by Hohenstein Research Institute.

Mains operation: 220-240 V; 50-60 Hz; 100 W

Product dimensions: approx. 62 x 42 cm

Product weight: approx. 580 g

3 year guarantee

Sales unit: 5 / Shipping carton: -

EAN no.: 4211125 21505 7

Item no.: 215.05

