




HK 72

beurer
wellbeing

HK 72 Akku Mobile Heat Belt

 german
engineering



4 illuminated temperature settings



Runtime of approx. 4 hours



Machine washable up to 30 °C



Up to 4 hours mobile heat independent from the mains - for maximum mobility

Powerful lithium-ion battery pack with LED displays for charge function and temperature settings

Economic shape for abdomen and back

Thin material, comfortable to wear under clothing

Cosy and skin friendly

Wide rubber strap with hook and loop fastener - comfortable to wear

4 illuminated temperature settings

Automatic switch-off after approx. 90 minutes

Rapid heat up

Temperature monitoring with safety switch-off

With charger (charging time approx. 2-3 hours)

Machine washable up to 30°C

4 W

Power supply:

Input: 100-240V, 50-60 Hz, 0.6 A

Output: 8.4V DC, 1.2 A

Product measurements (LxW): 115 x 18 cm

3 year guarantee

Sales unit: 3 / Shipping carton: -

EAN no.: 4211125 22591 9

Item no.: 225.91



225.91_HK72_2018-07-17_01_DS_EN_BEU Subject to errors and changes