




HK 123 XXL
Nordic

beurer
wellbeing

HK 123 XXL Nordic XXL heat pad



3
temperature
settings

3 illuminated temperature settings



cosy faux
fur fleece

Luxurious upper material -
fleecy soft and beautifully cuddly



machine
washable

Machine washable up to 30 °C



Extra large XXL-size: approx. 60 x 30 cm

Breathable, cosy and skin-friendly

Electronic temperature control

Automatic switch off after approx. 90 minutes

Rapid heat up

BSS (Beurer Safety System)

Detachable control unit

OEKO-TEX®

The textiles used for this device meet the stringent human ecological requirements of OEKO-TEX® Standard 100, as verified by Hohenstein Research Institute.

Mains operation:

220-240 V; 50-60 Hz; 100 W

Product dimensions (LxB): approx. 60 x 30 cm

Product weight: approx. 320 g

3 year guarantee

Sales unit: 5 / Shipping carton: -

EAN no.: 4211125 27416 0

Item no.: 274.16

