

MP 100

beurer
beauty



MP 100 Manicure/pedicure station

 german
engineering



Incl. 10 professional attachments
made from sapphire and felt



Speed is continuously adjustable,
clockwise/anti-clockwise rotation



Powerful professional motor –
max. power up to 20,000 rpm



For professional nail and foot care

Max. power up to 20,000 rpm

Steplessly adjustable

Runs in clockwise/counter-clockwise direction

10 high-quality sapphire and felt attachments

30 one-way attachments

Protective cap for nail dust

With non-slip handle surface

Safety lock

Suitable for diabetics

Easy to use

Also suitable for artificial nails

Mains operation

Includes integrated adapter storage

Speed: 5.000 rpm - 20.000 rpm \pm 20%

Mains operation:

Input: 100–240 V~, 50/60 Hz, 0.8 A

Output: 15.0 V— 2.0 A, 30.0 W

Product measurements (LxWxH): 22 x 23.5 x 11 cm

Product weight: approx. 1180 g

3 year guarantee

Sales unit: 3 / Shipping carton: -

EAN no.: 4211125 57068 2

Item no.: 570.68