



HT 15 Lice Comb


 german
 engineering


Simple combing frees hair from head lice and nits



No use of chemical agents

Acoustic signal stops by contact with a head lice or nit

The head lice and nits are removed by an electrical current upon contact

Including cleaning comb and protection cap

Medical Device

Incl. 1x 1.5V battery LR6

Product measures (LxWxH): 158 x 37 x 50 mm

Product weight: approx. 56 g (without battery)

3 year guarantee

Sales unit: 4 / Shipping carton: 4x4

EAN no.: 4211125 59168 7

Item no.: 591.68

