



HT 15 Lice Comb


**german
engineering**

**keep
lice free**

Simple combing frees hair from head lice and nits


**without
chemical
substances**

No use of chemical agents

Acoustic signal changes by contact with a head lice or nit

The head lice and nits are removed by an electrical current upon contact

Including cleaning comb and protection cap

Medical Device

Suitable for dogs and cats

Incl. 1x 1.5V battery LR6

Product measures (LxWxH): 158 x 37 x 50 mm

Product weight: approx. 56 g (without battery)

3 year guarantee

PZN: 16822224

Sales unit: 4 / Shipping carton: 4x4

EAN no.: 4211125 59334 6

Item no.: 593.34

