



IH 21 Nebuliser


 german
 engineering


Nebulisation with compressed
 air technology



For the treatment of the upper and lower
 respiratory tract

Nebulisation of liquid medicine using compressed air
 technology

For use in colds, asthma, respiratory disease, etc.

High proportion of respirable particles

High nebulisation capacity (approx. 0.3 ml/min)
 - short inhalation time

Accessories:

Nebuliser bottle, mouthpiece, nosepiece, adult mask, child mask,
 compressed air tube (extra long), filter

With practical storage compartment for the accessories

Can be disinfected

Medical product

Operating pressure: approx. 0.8 - 1.45 bar

230 V~ / 50 Hz / 230 VA

Product measurements: 30 x 18 x 10 cm

Product weight: approx. 1.65 kg

5 year guarantee

Sales unit: 2 / Shipping carton: 2 x 2

EAN no.: 4211125 60112 6

Item no.: 601.12