



IH 58 Nebuliser


 german
 engineering


Nebulisation of liquid medicine using compressed air technology



Very quiet, small and light

For the treatment of the upper and lower respiratory tract

For use in colds, asthma, respiratory disease, etc.

- for fast and effective inhalation

High proportion of respirable particles

High nebulisation capacity (approx. 0.25 ml/min.)

- short inhalation time

Accessories:

Nebuliser bottle, mouthpiece, nosepiece, adult mask, child mask, compressed air tube, mains adapter with micro-USB cable, filter, storage bag and carrying bag

Can be disinfected

Operating pressure: approx. 0.25-0.5 bar

Medical device

Mains operation:

Input: 100-240 V~; 50-60 Hz; 0.5 A

Output: 5 V $\overline{\text{---}}$ 2 A

Product measurements (LxWxH): 110 x 62 x 47 mm

Product weight: approx. 176 g

5 year guarantee

Sales unit: 3 / Shipping carton: 4 x 3

EAN no.: 4211125 60202 4

Item no.: 602.02

