



TL 20

**beurer**  
wellbeing

## TL20 Daylight therapy lamp

 **german  
engineering**



Simulation of daylight:  
Light intensity of approx. 10,000 lux (distance 10 cm\*)



Energy-saving light source



Medical device



For use in case of light deprivation symptoms in the winter months:  
e.g. low mood, lack of energy and drive

Compact size through LED technology

Small and handy - Ideal for the office or on the way

Exceptionally bright and even illumination

Light intensity of approx. 10,000 lux (distance 10 cm\*)

Flicker-free

UV-free

Energy-saving light source

Convenient one-button operation

Versatile stand

- Horizontal or vertical applicable
- Two tilt angles each way

Includes mains adapter, cable length approx. 175 cm

Certified medical device

Colour temperature: 6,000 K +/- 300 K

Illumination surface: approx. 16,5 x 8,5 cm

Product measurements: approx. 21 x 13 x 2 cm

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 608033

Item no.: 608.03



* Lux	Abstand / Distance / Distance Distanza / Distanza
5.000	ca./approx./env./aprox./ca. 20 cm
2.500	ca./approx./env./aprox./ca. 30 cm

608.03\_2018-03-15\_01\_EN Irrtum und Änderungen vorbehalten