



IL 35 Infrared lamp


**german
engineering**
150
watts

Output: 150 watts

medical
device

Certified medical device

auto
off3-level electronic timer with automatic
switch-off function

Soothing heat

For colds and muscle strains

Pressed glass bulb

3-level electronic timer with LED display

Switch-off function once set timer level has elapsed

Exclusive design

Adjustable shade: 5 angles

150 watts

Medical device

Product measurements: approx. 14.5 x 17 x 26 cm

Product weight: approx. 800 g

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 61612 0

Item no.: 616.12

