



## FM 150 Pro Vein Trainer



Air pressure massage of femur and lower legs with inflating and deflating air cushion



Continuously adjustable massage intensity



Improves the venous blood circulation

 **german engineering**



Revitalising pressure massage

Promotes blood circulation and relieves heavy and tired legs

Loosens tension

Air pressure massage with inflating and deflating air cushion

Continuously adjustable massage intensity

With 4 cuffs: 2 femur and 2 lower leg cuffs

Alternate inflation of the leg cuffs

Incl. Safety STOP button

With deflation function - for emptying the air cushions quickly

Incl. timer function 10/20/30 min.

Easy operation with manual control switch

Versatile application - whether you're at home or on the go

Individual adjustment of the leg cuff using hook-and-loop fasteners

Storage pouch included

Medical device

Battery or mains adaptor driven

Mains operation:

Input: 100-240V~; 50/60Hz; 200 mA

Output: 6V  $\overline{\overline{\overline{\quad}}}$  500 mA; 3 W

Incl. 4 x 1.5V battery AA

Product measurements (LxWxH): Femur cuff: 85 x 23,5 x 0,5 cm  
 Lower leg cuff: 61 x 67 x 0,5 cm

Product weight: approx. 1,7 kg

3 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 64327 0

Item no.: 643.27