



EM 26 Shoulder TENS

PAIN
 free


**german
 engineering**


TENS technology
 Transcutaneous electrical nerve stimulation (TENS)
 for drug free pain relief



water contact electrodes
 2x
 No contact gel, no replacement
 (electrodes with water contact)



4 training programs
 4 training programs



For shoulder pain relief

2 wear free water contact electrodes

Universal cuff for the left & right shoulder

Flexible ergonomic shape

Diagonal upper body circumference: up to approx. 120 cm

Controller with handy belt clip

Countdown timer

Safety switch-off

Battery changing indicator

Medical product

The device must not be used by users of cardiac pacemakers, epileptics, pregnant woman, people with cardiac arrhythmias

Output voltage max. 50V p-p (500 Ohm load)

Output current: max. 100 mA p-p (500 Ohm load)

Output frequency: 2-110 Hz

Pulse width: 60-220 µs per phase

Pulseform: symmetrical, biphasic square pulses

Incl. 3 x 1.5V batteries AAA

Product measurements (LxWxH): 8.6 x 8.3 x 3.6 cm

Product weight: approx. 90 g (incl. batteries)

Device, cuff and cable: 165 g (incl. batteries)

5 year guarantee

Sales unit: 3 / Shipping carton: 6 x 3

EAN no.: 4211125 64661 5

Item no.: 646.61

