




EM 27 Ankle TENS

PAIN
 free


**german
engineering**

**TENS
technology**

Transcutaneous electrical nerve stimulation (TENS)
for drug free pain relief


**water
contact
electrodes**
2x

No contact gel, no replacement
(electrodes with water contact)


**4 training
programs**

4 training programs



For pain relief on ankle

2 wear free water contact electrodes

Ankle cuff (adjustable with hook and loop fastener)

Flexible ergonomic shape

Circumference approx. 18-31 cm

Controller with handy belt clip

Countdown timer

Safety switch-off

Battery changing indicator

Medical product

The device must not be used by users of cardiac pacemakers,
epileptics, pregnant woman, people with cardiac arrhythmias

Output voltage max. 50V p-p (500 Ohm load)

Output current: max. 100 mA p-p (500 Ohm load)

Output frequency: 2-110 Hz

Pulse width: 60-220 μ s per phase

Pulseform: symmetrical, biphasic square pulses

Incl. 3 x 1.5V batteries AAA

Product measurements (LxWxH): 8.6 x 8.3 x 3.6 cm

Product weight: approx. 90 g (incl. batteries)

Cuff: approx. 60 g

5 year guarantee

Sales unit: 3 / Shipping carton: 6 x 3

EAN no.: 4211125 64707 0

Item no.: 647.07

