




FM 250
Vital Legs

beurer

FM 250 Vital Legs EMS circulation booster

 german
engineering



Electrostimulation device for relieving swelling and pain



Additional step function for mobilising the feet



Promotes circulation and relieves muscle tension

6 electrodes (2 electrode pads in the foot area and 1 electrode in the cuff)

2 separately adjustable channels

1 program with 15 pulse wave variations

Selectable intensity levels of 1 to 99

Incl. timer function (5 to 90 minutes)

Incl. cuffs for additional stimulation of the lower leg area

2 level locking function of the foot surface

Incl. hand switch for convenient use

Medical device

Incl. 4 x 1.5 V batteries AA and 2 x 1.5V batteries AAA

Mains operation:

Input: 100-240 V~; 50/60 Hz; 0.2 A

Output: DC 6.0 V === 300 mA

Product measurements (LxWxH): 37.2 x 28.7 x 11.1 cm

Product weight: approx. 1.5 kg (without batteries)

3 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 64937 1

Item no.: 649.37



649.37_FM250_VitalLegs_2020-10-23_02_DS_EN_BEU Subject to errors and changes

PhysioLine