



BC 21

beurer
medical

BC 21 wrist blood pressure monitor

 **german
engineering**



Voice output in five different languages
(voice output disengageable):
(A): EN, ES, PL, AR, FA;
(B): DE, EN, FR, IT, TR



Risk indicator:
Classification of measurements
with coloured scale



Fully automatic blood pressure and pulse measurement on the wrist
With cuff for wrist circumferences from 14.0 - 19.5 cm

Arrhythmia detection: Warning in case of possible heart rhythm disturbances

Two user memories for 60 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Clear, easy-to-read display

Adjustable volume

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 2 x 1.5V AAA batteries and storage case

Medical device

Product measurements (LxWxH): 87 x 67 x 31 mm

Product weight: approx. 105 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: (A): 4211125 65040 7

(B): 4211125 65046 9

Item no.: (A): 650.40

(B): 650.46