



BC 40

beurer  
medical

## BC 40 wrist blood pressure monitor



Fast and convenient measurement, even during pumping up



Risk indicator: Classification of measurements with coloured scale



Fully automatic blood pressure and pulse measurement on the wrist

With cuff for wrist circumferences from 12.5 - 21.5 cm

Arrhythmia detection: Warning in case of possible heart rhythm disturbances

User memory for 60 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Clear, easy-to-read display

Date and time

Automatic switch-off

Signal in case of application error

Battery changing indicator

Incl. 2 x 1.5V AAA batteries and practical storage case

Medical device

Product measurements (LxWxH): 76 x 70 x 23 mm

Product weight: approx. 105 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65061 2

Item no.: 650.61