



BM 26

beurer
medical

BM 26 – Blood pressure monitor

 german
engineering

risk
indicator



Risk indicator:
Classification of measurements with coloured scale

irregular
heart rhythm



Arrhythmia detection: Warning in case of possible
heart rhythm disturbance



Fully automatic blood pressure and pulse measurement on the upper arm

Four user memories for 30 measured values each

With cuff for upper arm circumferences from 22-35 cm

Average value calculation of the last 3 measurements

Clear, easy-to-read display

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 4 x 1,5V AA batteries and storage pouch

Medical device

Product measurements (LxWxH): 155 x 110 x 70 mm

Product weight: approx. 395 g (without batteries, with cuff)

XL cuff (30-42 cm) available separately: Art. No. 162.973

5 year guarantee

Sales unit: 4 / Shipping carton: 6 x 4

EAN no.: 4211125 65228 9

Item no.: 652.28