



BM 44


beurer
medical

BM 44 – Blood pressure monitor

 german
engineering



Easy to use via one single button



Blue illuminated,
easy-to-read XL display



Fully automatic blood pressure and pulse measurement on the upper arm

Risk indicator: Classification of measurements with coloured scale

Illuminated START/STOP button

Arrhythmia detection: Warning in case of possible heart rhythm disturbance

With cuff for upper arm circumferences from 22 - 30 cm

Can be used with mains part (available separately)

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Battery changing indicator

Incl. 4 x 1.5V AA batteries and practical storage pouch

Medical device

Product measurements (LxWxH): 94 x 122 x 53 mm

Product weight: approx. 442 g (without batteries, with cuff)

XL cuff (30 - 42 cm) available separately: Art. No. 162.795

Mains part available separately: Art. No. 071.95

5 year guarantee

Sales unit: 4 / Shipping carton: 3 x 4

EAN no.: 4 211125 65501 3

Item no.: 655.01