



BM 54

beurer  
medical

## BM 54 Bluetooth® Blood pressure monitor

 german  
engineering



With cuff for upper arm circumferences  
from 22 - 44 cm



Classification of measurements  
with coloured scale



Fully automatic blood pressure and pulse measurement on the upper arm

Arrhythmia detection: Warning in case of possible  
heart rhythm disturbances

Two user memories for 60 measured values each

Average value of all stored measurements and morning and evening blood  
pressure for the last 7 days

Easy-to-read XL display

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 4 x 1,5V AAA batteries and storage pouch

Medical device

Product measurements (LxWxH): 139 x 94 x 48 mm

Product weight: approx. 349 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65512 9

Item no.: 655.12



## CONNECTION



Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.



Syncs with  
 • iOS ≥ 8.0  
 • Android™ ≥ 4.4  
 • Bluetooth® ≥ 4.0



List of supported devices

\* tested online application "beurer HealthManager" Cloud  
 Can only be used with its full functionality with the „beurer HealthManager“ app.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.