



BC 87

beurer
medical

BC 87 Bluetooth® Wrist blood pressure monitor



Patented resting indicator* checks the resting state during the measurement



Inflation Technology: Fast and convenient measurement, even during pumping up



Accurate measurement results due to correct measurement position at heart level



Fully automatic blood pressure and pulse measurement on the wrist

With cuff for wrist circumferences from 13,5 - 21,5 cm

Clear, easy-to-read black XL display

Extra slim design – the ideal companion, even on the go

Risk indicator: Classification of measurements with coloured scale

Irregular heartbeat detection

Two user memories for 120 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Battery changing indicator

Incl. 2 x 1,5V AAA batteries and practical storage case

Medical device

Product measurements (LxWxH): 72 x 96 x 71 mm

Product weight: approx. 119 g (without batteries, with cuff)

5 year guarantee


Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4 211125 657031

Item no.: 657.03

*Patent EP1673 009

657.03_BC87_2021-08-11_01_DS_EN_BEU | Subject to errors and changes




beurer Grafik

Released without change
 Released after change
 Change and resubmit

25.08.21 jaan
 Date Signature



VERBINDUNG



Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.



Syncs with

- iOS ≥ 12.0
- Android™ ≥ 8.0
- Bluetooth® ≥ 4.0



List of supported devices

657.03_BC87_2021-08-11_01_DS_EN_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.