



BC 85

beurer  
medical

## BC 85 Bluetooth® Wrist blood pressure monitor

 german  
engineering



For correct positioning of the blood pressure monitor at heart level



Fast and convenient measurement, even during pumping up



Large, perfectly readable XL display



Fully automatic blood pressure and pulse measurement on the wrist

With cuff for wrist circumferences from 13.5 - 23.0 cm

Innovative networking between smartphone and blood pressure monitor via Bluetooth® low energy technology

USB port

Risk indicator: Classification of measurements with coloured scale

Arrhythmia detection: Warning in case of possible heart rhythm disturbances

Extra slim design - the ideal companion, even on the go

Two user memories for 60 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 2 x 1.5V AAA batteries, USB cable and practical storage box

Medical device

Product measurements (LxWxH): 97 x 68 x 22 mm

Product weight: approx. 108 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65706 2

Item no.: 657.06



**Positioning indicator**

- For correct positioning of the blood pressure monitor at heart level



**Risk indicator**

- Classification of measurements with coloured scale




**USB port**

- USB interface for innovative networking between the blood pressure monitor and "beurer HealthManager" PC software



CONNECTION



\* tested online application "beurer HealthManager" Cloud  
Can only be used with its full functionality with the „beurer HealthManager“ app.

657.06\_BC85\_2019-01-25\_01\_DS\_EN\_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.