



BM 85

beurer
medical

BM 85 Bluetooth® Blood pressure monitor

german
engineering



Resting indicator* for exact measurement result



With cuff for upper arm circumferences
from 22 - 42 cm



Fully automatic blood pressure and pulse
measurement on the upper arm

Inflation Technology: Fast and convenient measurement,
even during pumping up

LED risk indicator: Coloured LED scale for classifying
measurements

Arrhythmia detection: Warning in case of possible heart
rhythm disturbances

2 x 60 memory spaces

Average value of all stored measurements and morning and
evening blood pressure for the last 7 days

Easy to use via illuminated sensor touch buttons

White illuminated XL display

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. USB power mains part, cuff holder and storage pouch

With lithium-ion battery

Medical device

Mains operation:

Input: 100-240 V~; 50-60 Hz; 60-120 mA

Product measurements (LxWxH): 180 x 99 x 40 mm

Product weight: approx. 456 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 2 x 4

EAN no.: 4211125 65803 8

Item no.: 658.03

*Patent EP 1 673 009

658.03_BM85_2019-02-04_02_DS_EN_BEU Subject to errors and changes



LED risk indicator

• Coloured LED scale for classifying measurements



Sensor touch

• Easy to use via illuminated sensor touch buttons




USB port

• USB interface for innovative networking between the blood pressure monitor and "beurer HealthManager" PC software and charging the lithium-ion battery



CONNECTION



Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.



Syncs with
• iOS ≥ 10.0
• Android™ ≥ 4.4
• Bluetooth® ≥ 5.0



List of supported devices

658.03_BM85_2019-02-04_02_DS_EN_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

* tested online application "beurer HealthManager" Cloud
Can only be used with its full functionality with the „beurer HealthManager“ app.