



BC 32 wrist blood pressure monitor


 german
 engineering


Risk indicator: Classification of measurements with coloured scale



Arrhythmia detection:
Warning in case of possible heart rhythm disturbances



Fully automatic blood pressure and pulse measurement on the wrist

With cuff for wrist circumferences from 13.5 - 19.5 cm

Two user memories for 60 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Clear, easy-to-read display

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 2 x 1.5V AAA batteries and practical storage case

Medical device

Product measurements (LxWxH): 72 x 69 x 29.2 mm

Product weight: approx. 101 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 4211125 65909 7

Item no.: 659.09