



BC 58

beurer
medical

BC 58 Wrist blood pressure monitor

german
engineering

beurer
CONNECT



sensor
touch

Easy to use via sensor touch buttons



xl
display

Clear, easy-to-read XL black display



usb
port

For innovative networking between
the blood pressure monitor and
"beurer HealthManager" PC software



Fully automatic blood pressure and pulse measurement on the wrist

With cuff for wrist circumferences from 14 - 19.5 cm

Risk-indicator: Classification of measurements with coloured scale

Arrhythmia detection: Warning in case of possible heart rhythm disturbances

Two user memories for 60 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Battery changing indicator

Incl. USB cable, 2x 1.5 AAA batteries and practical storage case

Medical device

Product measurements (LxWxH): 90 x 68 x 30 mm

Product weight: approx. 158 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65916 5

Item no.: 659.16



Sensor touch

- Easy to use via sensor touch buttons



USB port

- USB interface for innovative networking between the blood pressure monitor and „beurer HealthManager“ PC software



CONNECTION



SOFTWARE

For innovative networking the blood pressure monitor and „beurer HealthManager“



* tested online application „beurer HealthManager“ Cloud
Can only be used with its full functionality with the „beurer HealthManager“ app.