




BC 58

beurer
medical

BC 58 Wrist blood pressure monitor



beurer
CONNECT



Easy to use via sensor touch buttons



Clear, easy-to-read XL black display



For innovative networking between the blood pressure monitor and free web application incl. USB uploader "beurer HealthManager Pro" / PC software "beurer HealthManager"



Fully automatic blood pressure and pulse measurement on the wrist

With cuff for wrist circumferences from 14 - 19.5 cm

Risk-indicator: with integrated scale for grading the measurement results in the display

Arrhythmia detection: Warning in case of possible heart rhythm disturbances

Two user memories for 60 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Battery changing indicator

Incl. USB cable, 2x 1.5 AAA batteries and practical storage case

Medical device

Product measurements (LxWxH): 90 x 68 x 30 mm

Product weight: approx. 158 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65916 5

Item no.: 659.16



Sensor touch

- Easy to use via sensor touch buttons




USB port

- USB interface for innovative networking between the blood pressure monitor and free web application incl. USB uploader "beurer HealthManager Pro" / PC software "beurer HealthManager"



CONNECTION



SOFTWARE

For innovative networking between the blood pressure monitor and free web application incl. USB uploader "beurer HealthManager Pro" / PC software "beurer HealthManager"



HealthManager Pro



HealthManager



*

* tested online application „beurer HealthManager“ Cloud