




EM 38 Back belt with TENS device

PAIN
 free


**german
 engineering**


Transcutaneous electrical nerve stimulation (TENS)



4 non-wearing contact electrodes of
 conductive carbon material



Easy to use

With velcro fastening

Stomach belt for waist sizes from 30 to max. 55 in
 (75 to max. 140 cm), extension belt included

Flexible ergonomic shape

Nonwearing contact electrodes made of
 conductive carbon material

Electrodes with water contact (no contact gel, no replacement)

LCD display with figurative symbols

4 training programs (25-30 min)

Intensity adjustable from 0 up to 20

Countdown timer with remaining time

Safety switch-off / pause function / automatic switch-off / key lock

Removable control element

Belt hand washable

Battery changing indicator

The device must not be used by users of cardiac pacemakers,
 epileptics, pregnant woman, people with cardiac arrhythmias

Output voltage max. 50V p-p (500 Ohm load)

Output current: max. 100 mA p-p (500 Ohm load)

Output frequency: 2-110 Hz

Pulse width: 60-220 μ s per phase

Pulseform: symmetrical, biphasic square pulses

Batteries included

Medical product

5 year guarantee

Sales unit: 3 / Shipping carton: 6 x 3

EAN no.: 4211125 664008

Item no.: 664.00

