



SL 40 Snore belt


**german
engineering**
anti-snoring


To prevent snoring that results from lying on your back

**EMS
technology**


You are immediately prompted to stop lying on your back through electrostimulation on the back

**adjustable
intensity**


Adjustable intensity



For the prevention of snoring that results from lying on your back through a position sensor that recognises when you are on your back

You are immediately prompted to stop lying on your back through electrostimulation on the back

To make you change your sleeping position

For light to heavy snoring

Easy to use

Attachment around the chest and back

Comfortable shape

To be worn unobtrusively under clothing

Chest circumference approx. 80 - 105 cm

Adjustable intensity

Varying impulses in order to prevent you getting used to them

Adjustable length

Battery operated

- Charging time approx. 1 hour

- Runtime approx. 2 - 5 nights

Removable control unit on the belt fastener

Strap can be hand-washed

2 wear-free contact electrodes made of conductive carbon material

No contact gel or replacement electrodes required (water contact electrodes)

Safety switch-off

Incl. extension belt (chest circumference approx. 140 cm)

Incl. micro USB charging cable

Product measurements: approx. 117 x 6.5 x 1.5 cm

Product weight: approx. 175 g

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 67749 7

Item no.: 677.49

Control unit

Extension belt
