



BF 220 Diagnostic bathroom scale


 german
 engineering


Large platform made from
8 mm safety glass



Blue illuminated XL-display
Digit size: 40 mm



With brushed stainless steel electrodes
Determination of body weight, body fat,
body water and muscle ratio

Graduation:

- 0.1% for body fat, body water and muscle ratio
- 100 g for body weight

Change over to kg/lb/st

Capacity: 180 kg

5 activity levels

Quick-Start-Technology

Overload indication

Automatic switch-off

Incl. 2 x 3V battery CR2032

Product measurements (LxWxH): 32 x 32 x 2.5 cm

Product weight: approx. 2400 g

5 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 74813 5

Item no.: 748.13