



# BF 700 Diagnostic bathroom scale

 **german engineering**



Illuminated LCD display



Calculates body fat, body water, muscle percentage, bone mass, BMI, BMR/AMR (kcal)



With app + software for downloading: analysis records and graphical evaluations of body values



- Innovative networking between smartphone and scale
- With *Bluetooth*<sup>®</sup> low energy technology - keep track of your measurements and body data at home and on the move
- inclusive "beurer HealthManager" and "beurer BodyShape" App for downloading
- Automatically data transfer within the *Bluetooth*<sup>®</sup> range if the app is open
- Black LCD display with white illuminated digits
- Digit size: 23 mm
- Displays 3 initials of the user
- Electrodes made of brushed stainless steel
- Display of body weight, body fat, body water, muscle ratio, bone mass, AMR/BMR and BMI
- Graduation:
  - 0,1% body fat, body water, muscle ratio
  - 100 g for body weight, bone mass
- 5 levels of activity
- Capacity: 180 kg
- 8 user memories with automatic user recognition
- Automatic switch off
- Quickstart
- Change over to kg/lb/st
- Overload indication
- 3 x 1.5V AAA batteries included
- Product measurements (LxWxH): 30 x 30 x 2.3 cm
- Product weight: approx. 2040 g
- 5 year guarantee
- Sales unit: 3 / Shipping carton: -
- EAN no.: 4211125 74834 0
- Item no.: 748.34

**body analysis**



**weight**



**body fat**



**body water**



**muscle percentage**



**bone mass**



**BMI**



**BMR/AMR (kcal)**






 **german  
engineering**



## CONNECTION

**beurer connect**



Innovative networking between smartphone and scale - keep track of your measurements and body data at home and on the move. Easily synchronised with your smartphone via *Bluetooth®*

---

works with



Download on the **App Store**

GET IT ON **Google Play**

System requirements:  
 iOS ≥ 10.0  
 Android™ ≥ 5.0  
 Bluetooth® ≥ 4.0



List of supported smartphones



\*\*\* tested online application "beurer HealthManager" Cloud

748.34\_BF700\_2019-07-18\_01\_DS\_EN\_BEU Subject to errors and changes

The *Bluetooth®* word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

Android is a trademark of Google LLC.