




BF 950

beurer
wellbeing

BF 950 Diagnostic bathroom scale


german
engineering



Colour LED response for achieving targets



Pregnancy mode for controlled weight gain during pregnancy (without BIA*)



ITO** coated weighing surface made of safety glass
Automatic data transfer within Bluetooth® range after opening the app

Determination of body weight, body fat, body water, muscle percentage, bone mass, AMR/BMR and BMI

XL display in modern dot matrix look

Digit size: 23 mm

Time/temperature display

Display of user initials

Change over to kg/lb/st

Capacity: 180 kg

Graduation:

- 0,1% for body fat, body water and muscle ratio

- 100 g for body weight and bone mass

5 activity levels

8 user memories with automatic user recognition

Tare function

Overload indication

Quick-Start-Technology

Automatic switch-off

Incl. 4 x 1.5V batteries AAA

Product measurements (LxWxH): 32 x 32 x 2.4 cm

Product weight: approx. 2.1 kg

5 year guarantee

Sales unit: 3 / Shipping carton: -

EAN no.:

Black: 4211125 74912 5

White: 4211125 74911 8

Item no.:

Black: 749.12

White: 749.11

* Bioelectrical impedance analysis
** Indium Tin Oxide




BF 950

beurer
wellbeing

german
engineering



CONNECTION



Innovative networking between smartphone and scale - keep track of your measurements and body data at home and on the move. Easily synchronised with your smartphone via *Bluetooth*®



Syncs with
 • iOS ≥ 8.0
 • Android™ ≥ 4.4
 • Bluetooth® ≥ 4.0



List of supported devices

748.11-749.12_BF950_2018-06-15_01_DS_EN_BEU Subject to errors and changes

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

The *Bluetooth*® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.