




FT 09/1 Thermometer blue



waterproof

100% waterproof



auto off

Automatic switch off

including
protective
cover

With protective cover



Digital thermometer

Large Display

3-digit display

Audible signal when measuring stops

1 memory place

Measurement accuracy $\pm 0.1^{\circ}\text{C}$ (between 35.5 and 42°C)Measurement accuracy $\pm 0.2^{\circ}\text{C}$ (between $32 - 35.4^{\circ}\text{C}$ and $42.1 - 42.9^{\circ}\text{C}$)

100% waterproof

Can be disinfected

Without glass, mercury free

Incl. 1 x 1,5V battery LR41

Battery changing indicator

Product measurements: 129 x 20 x 17 mm

Product weight: approx. 10 g (including battery)

Medical device

5 year guarantee

Sales unit: 12 / Shipping carton: 12 x 12

EAN no.: 4211125 79116 2

Item no.: 791.16