



FT 58 Ear thermometer


 german
 engineering


Hygienic and safe measurement



With interchangeable protective cap
 Protective cap replacement set
 available



Measurement in seconds



Infrared measurement technology

Temperature alarm at 38°C

Warning sound at 37.5°C

Indication of °C and °F

Object temperature measurement

Room temperature display

10 memory spaces

Acoustic signal when measurement is complete

Laughing Smiley = everything is OK

Sad Smiley = $\geq 38^{\circ}\text{C}$ temperature alarm

Automatic switch off

With low battery indicator

With storage box

Incl. 10 protective caps

Incl. 2 x 1.5V batteries AAA

Product measurements (LxBxH): 155 x 38 x 47 mm

Product weight: approx. 63 g (without batteries)

Medical device

5 year guarantee

Sales unit: 4 / Shipping carton: 12 x 4

EAN no.: 4211125 79533 7

Item no.: 795.33