



 german|engineering

## BF 180 Diagnostic bathroom scale

blue  
illuminated  
display



Blue illuminated display

large  
display  
34 mm



With extra large display

Digit size: 34 mm



Determination of body weight, body fat, body water,  
muscle ratio, bone mass and calory consumption (BMR)

Calculation of BMI for simple classification of your  
body weight

With brushed stainless steel electrodes

Platform made of safety glass

Change over to kg/lb/st

Capacity: 180 kg

Graduation:

- 0.1% for body fat, body water and muscle ratio

- 100 g for body weight and bone mass

5 activity levels

10 user memories

Overload indication

Quick-Start-Technology

Automatic switch-off

Incl. 2 x 3V batteries CR2032

Product measurements (LxWxH): 30 x 30 x 1.9 cm

Product weight: approx. 1.55 kg

5 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 74909 5

Item no.: 749.09