



\*



BF 850

beurer  
wellbeing

german|engineering

## BF 850 Diagnostic bathroom scale



Innovative networking between smartphone and scale  
With free "beurer HealthManager" app + software for downloading



With *Bluetooth*® low energy technology  
- Keep track of your measurements and body data  
at home and on the move



Particularly large, modern framed platform



Automatic data transfer within *Bluetooth*® range after opening the app

Determination of body weight, body fat, body water, muscle percentage, bone mass, AMR/BMR and BMI

Large 2 lined LCD display with blue digit illumination

Digit size: 22 mm

Display of user initials

Incl. "beurer HealthManager" App for downloading

With brushed stainless steel electrodes

Platform made of safety glass

Change over to kg/lb/st

Capacity: 180 kg

Graduation:

- 0,1% for body fat, body water and muscle ratio

- 100 g for body weight and bone mass

5 activity levels

8 user memories with automatic user recognition

Overload indication

Quick-Start-Technology

Automatic switch-off

Incl. 3 x 1.5V batteries AAA

Product measurements (LxWxH): 32.5 x 32.5 x 2.4 cm

Product weight: approx. 2 kg

5 year guarantee

Sales unit: 3 / Shipping carton: -

EAN no.:

Black: 4211125 74821 0

White: 4211125 74824 1

Item no.:

Black: 748.21

White: 748.24

\* tested online application "beurer HealthManager" Cloud

0717 Subject to errors and changes



\*




BF 850

beurer  
wellbeing

german|engineering



## CONNECTION



Innovative networking between smartphone and scale - keep track of your measurements and body data at home and on the move. Easily synchronised with your smartphone via *Bluetooth®*



Syncs with

- iOS ≥ 8.0
- Android™ > 4.4
- Bluetooth® 4.0



List of supported devices