



germanengineering

BG 55 Glas-Diagnosewaage

Designer-Glasdiagnosewaage

XXL-Display (18 x 4 cm) ermöglicht simultanes Betrachten der Ergebnisse (Gewichtsanzeigen- Zifferngröße 28 mm)

Extra große Trittpläche aus Sicherheitsglas (ca. 32 x 32 cm)

Elektroden aus sandgestrahltem Edelstahl

Ermittlung von Körpergewicht, Körperfett, Körperwasser, Muskelanteil und Knochenmasse

Einteilung:

- 0,1% für Körperfett, Körperwasser und Muskelanteil
- 100 g für Körpergewicht und Knochenmasse

Interpretation von Körpergewicht, Körperfett, Körperwasser und Muskelanteil in feingliedriger Aufteilung

Trendanzeige (Messwert wird mit dem Durchschnitt der letzten 5 Messwerte verglichen)

BMR + AMR (Grundumsatz und Aktivitätsumsatz in Kcal)

Sportlermodus mit 5 Aktivitätsgraden

150 kg Tragkraft

10 Benutzerspeicher und jeweils 5 Messwertspeicher

Abschaltautomatik

Überlastungsanzeige

Tap-On Technik

Umstellung von kg/lb/st

Produktmaße: 35 x 32 x 3,4 cm

Produktgewicht: ca. 2200 g

Inklusive 3 x 1,5V AAA Batterien

5 Jahre Garantie

VE: 3 / Transportkarton:

EAN-Nr.: 4211125 75018 3

Art.-Nr.: 750.18



BG 55 Glass diagnostic scale

Design glass diagnostic scale

XXL display (18 x 4 cm): all results at one glance (digit size weight 28 mm)

Extra large platform of safety glass (approx. 32 x 32 cm)

Electrodes of sand blasted stainless steel

Display of body weight, body fat, body water, muscle ratio and bone mass

Division:

- 0.1% for body fat, body water and muscle ratio
- 100 g for body weight and bone mass

Interpretation of body weight, body fat, body water and muscle ratio in finest segmentation

Trend index (current result will be compared with the average of the last 5 results)

BMR + AMR (basic metabolic rate and active metabolic rate in kcal)

Athlete mode with 5 activity levels

150 kg capacity

10 user memories with 5 data memories each

Automatic switch off

Overload indication

Tap-on technology

Change over to kg/lb/st

Product measurements: 35 x 32 x 3.4 cm

Product weight: approx. 2200 g

3 x 1.5V AAA batteries included

5 year guarantee

Sales unit: 3 / Shipping carton:

EAN no.: 4211125 75018 3

Item no.: 750.18

