

Rückseite
Back side



EM 20 EMS-Pad Sixpack

EMS (Elektrische Muskel-Stimulation) - für Muskeltraining und -regeneration

Soft und flexibel, komfortabel und sicher

Punktgenaue Anwendung

Intensität einstellbar (15 Stufen)

Abschaltautomatik

Auswechselbare Batterie

Inkl. selbsthaftendem Gel-Film-Set

Inkl. 1 x 3V Batterie CR2032

Produktmaße: 170 x 170 x 13 mm

Produktgewicht: ca. 69 g (mit Batterie)

3 Jahre Garantie

Nachkaufset EM 20 EMS-Pad Sixpack

1 x Gel-Film-Set (selbsthaftend)

1 x 3V Batterie CR2032

VE: 10 / Transportkarton: 10

EAN-Nr.: 4211125 64716 2

Art.-Nr.: 647.16

VE: 10 / Transportkarton: 8 x 10

EAN-Nr.: 4211125 64711 7

Art.-Nr.: 647.11

EM 20 EMS Sixpack Pad

EMS (Electronic muscle stimulation) - for muscle training and regeneration

Soft and flexible, comfortable and safe

Targeted usage

Adjustable intensity (15 levels)

Automatic switch-off

Battery can be replaced

Incl. self-adhesive gel film set

Incl. 1 x 3V battery CR2032

Product measurements: 170 x 170 x 13 mm

Product weight: approx. 69 g (with battery)

3 year guarantee

Replacement set EM 20 EMS Sixpack Pad

1 x Gel-film-set (self-adhesive)

1 x 3V battery CR2032

VE: 10 / Shipping carton: 10

EAN no.: 4211125 64716 2

Item no.: 647.16

Sales unit: 10 / Shipping carton: 8 x 10

EAN no.: 4211125 64711 7

Item no.: 647.11