



EM 37 Bauchmuskel-Gürtel

Stimulation der zentralen und seitlichen Bauchmuskulatur

4 Elektroden

Einfache Anwendung

Flexibler Bauchgurt mit Klettverschluss

Universell für Bauchumfänge von 70 bis max. 140 cm verwendbar

Anschmiegsame ergonomische Form

Verschleißfreie Kontaktelektroden aus leitendem Carbonmaterial

Wasserkontakt-Elektroden (kein Kontaktgel oder Ersatzelektroden notwendig)

LCD Anzeige mit anschaulicher Symbolik

5 Trainingsprogramme (22-31 Minuten)

Intensität von 0 bis 40 einstellbar

Countdown-Timer mit Restlaufzeit

Trainingsspeicher-Funktion für fortlaufende Trainingskontrolle

Sicherheitsabschaltung/ Pause-Funktion/ Abschaltautomatik/ Tastensperre

Abnehmbares Bedienelement

Gurt handwaschbar

Batteriewechselanzeige

Die Anwendung ist verboten bei Trägern von Herzschrittmachern, Epileptikern, Schwangeren, Personen mit Herzrhythmusstörungen

Ausgangsspannung: max. 70V p-p bei 500 Ohm Ladung

Ausgangsstrom: max. 140 mA p-p bei 500 Ohm Ladung

Ausgangsfrequenz: 30-80 Hz

Pulsbreite: 200 µs pro Phase

Impulsform: symmetrische, zweiphasige Rechteckform

Inkl. 3 x 1,5V Batterien AAA

Produktmaße (LxBxH): 8,8 x 9,1 x 2,2 cm

3 Jahre Garantie

VE: 3 / Transportkarton: 6 x 3

EAN-Nr.: 4211125 64703 2

Art.-Nr.: 647.03

EM 37 Abdominal Muscle Belt

Stimulation of front and side abdominal muscles

4 electrodes

Easy to use

Flexible abdominal belt with hook and loop fastener

Flexible stomach belt for waist sizes from 28 to max. 55 in (70 to max. 140 cm)

Flexible ergonomic shape

Nonwearing contact electrodes made of conductive carbon material

Electrodes with water contact (no contact gel, no replacement)

LCD display with figurative symbols

5 training programs (22 to 31 min)

Intensity adjustable from 0 up to 40

Countdown timer with remaining time

Training memory function for continuous training control

Safety switch-off/ pause function/ automatic switch-off/ key lock

Removable control element

Belt hand washable

Battery changing indicator

The device must not be used by users of cardiac pacemakers, epileptics, pregnant woman, people with cardiac arrhythmias

Output voltage max. 70V p-p (500 Ohm load)

Output current: max. 140 mA p-p (500 Ohm load)

Output frequency: 30-80Hz

Pulse width: 200 µs per phase

Pulseform: symmetrical, dual-phase square pulses

Incl. 3 x 1.5V batteries AAA

Product measurements (LxWxH): 8.8 x 9.1 x 2.2 cm

3 year guarantee

Sales unit: 3 / Shipping carton: 6 x 3

EAN no.: 4211125 64703 2

Item no.: 647.03