



germanengineering

### EM 38 Elektrostimulationsgerät

Einfache Anwendung
Mit Klettverschluss
Universell für Bauchumfänge von 75 bis max. 140 cm verwendbar, inkl. Verlängerung
Anschmiegsame ergonomische Form
Verschleissfreie Kontaktelektroden aus leitendem Carbonmaterial
Wasserkontakt-Elektroden (kein Kontaktgel oder Ersatzelektroden notwendig)
LCD Anzeige mit anschaulicher Symbolik
4 Trainingsprogramme (25-30 Minuten)
Intensität von 0 bis 20 einstellbar
Countdown-Timer mit Restlaufzeit
Sicherheitsabschaltung/ Pause-Funktion/ Abschaltautomatik/ Tastensperre
Abnehmbares Bedienelement
Gurt handwaschbar
Batteriewechselanzeige
Die Anwendung ist verboten bei Trägern von Herzschritt-machern, Epileptikern, Schwangeren, Personen mit Herz-rhythmusstörungen
Ausgangsspannung: max. 50V p-p bei 500 Ohm Ladung
Ausgangsstrom: max. 100 mA p-p bei 500 Ohm Ladung
Ausgangsfrequenz: 2-110 Hz
Pulsbreite: 60-220 µs pro Phase
Impulsform: symmetrische, biphasische Rechteckform
Inklusive Batterien
Medizinprodukt
3 Jahre Garantie

VE: 3 / Transportkarton: 6 x 3

EAN-Nr.: 4211125 66400 8

Art.-Nr.: 664.00



### EM 38 Digital electronic stimulation unit

Easy to use
With velcro fastening
Stomach belt for waist sizes from 30 to max. 55 in (75 to max. 140 cm), extension belt included
Flexible ergonomic shape
Nonwearing contact electrodes made of conductive carbon material
Electrodes with water contact (no contact gel, no replacement)
LCD display with figurative symbols
4 training programs (25-30 min)
Intensity adjustable from 0 up to 20
Countdown timer with remaining time
Safety switch-off/ pause function/ automatic switch-off/ key lock
Removable control element
Belt hand washable
Battery changing indicator
The device must not be used by users of cardiac pacemakers, epileptics, pregnant woman, people with cardiac arrhythmias
Output voltage max. 50V p-p (500 Ohm load)
Output current: max. 100 mA p-p (500 Ohm load)
Output frequency: 2-110 Hz
Pulse width: 60-220 µs per phase
Pulseform: symmetrical, biphasic square pulses
Batteries included
Medical product
3 year guarantee

Sales unit: 3 / Shipping carton: 6 x 3

EAN no.: 4211125 664008

Item no.: 664.00

