

## GS 485 Glaswaage

Optimale Gewichts-Kontrolle via App + Software zum Download

Innovative Vernetzung zwischen Smartphone und Waage - Messwerte und Körperdaten zuhause und unterwegs immer im Blick

Übertragung per *Bluetooth*® Technologie

Komfortable Auswertung über „beurer HealthManager“ App und optionaler HealthManager Software

Elegant gerahmte Glaswaage

Wiegefläche aus Sicherheitsglas

LCD-Display

Zifferngröße: 36 mm

BMR + AMR (Grundumsatz und Aktivitätsumsatz in Kcal)

Sportlermodus mit 5 Aktivitätsgraden

180 kg Tragkraft

8 Benutzerspeicher mit automatischer Benutzererkennung

Abschaltautomatik

Quickstart

Umstellung von kg/lb/st

Überlastungsanzeige

Inklusive 1 x 3V CR2032 Batterie

Produktmaße (LxBxH): 35,5 x 32,0 x 2,2 mm

Produktgewicht: ca. 1,8 kg

5 Jahre Garantie

VE: 3 / Transportkarton: -

EAN-Nr.: 4211125 75612 3

Art.-Nr.: 756.12



## GS 485 Glass scale

Optimum weigh control via app + software for downloading

Innovative networking between smartphone and scale - keep track of your measurements and body data at home and on the move

Transfer via *Bluetooth*® technology

Convenient evaluation using the “beurer HealthManager“ app and the optional HealthManager software

Elegant framed glass scale

Weighing surface made of safety glass

LCD display

Digit size: 36 mm

BMR + AMR (basic metabolic rate and active metabolic rate in kcal)

Athlete mode with 5 activity levels

180 kg capacity

8 user memories with automatic user recognition

Automatic switch off

Quickstart

Change over to kg/lb/st

Overload indication

1 x 3V CR2032 battery included

Product measurements (LxWxH): 35.5 x 32.0 x 2.2 mm

Product weight: approx. 1.8 g

5 year guarantee

Sales unit: 3 / Shipping carton: -

EAN no.: 4211125 75612 3

Item no.: 756.12

