



germanengineering

## HK 35 Heizkissen

- \_\_\_\_\_ Schnellheizung
- \_\_\_\_\_ Abschaltautomatik nach ca. 90 Minuten
- \_\_\_\_\_ Elektronische Temperaturregelung
- \_\_\_\_\_ 3 Temperaturstufen
- \_\_\_\_\_ Beleuchtete Schaltstufen
- \_\_\_\_\_ Abnehmbare Zuleitung
- \_\_\_\_\_ BSS (Beurer Sicherheits System)
- \_\_\_\_\_ Abnehmbarer Baumwollbezug, waschbar
- \_\_\_\_\_ Öko-Tex 100
- \_\_\_\_\_ Die bei diesem Gerät eingesetzten Textilien erfüllen die hohen humanökologischen Anforderungen des Öko-Tex Standards 100, wie durch das Forschungsinstitut Hohenstein nachgewiesen.
- \_\_\_\_\_ 100 Watt
- \_\_\_\_\_ Produktmaße: 40 x 30 cm
- \_\_\_\_\_ Produktgewicht: ca. 460 g

3 Jahre Garantie

VE: 5 / Transportkarton:  
 EAN-Nr.: 4211125 26001 9  
 Art.-Nr.: 260.01



## HK 35 Heating pad

- \_\_\_\_\_ Rapid heat up
- \_\_\_\_\_ Automatic switch off after approx. 90 minutes
- \_\_\_\_\_ Electronic temperature control
- \_\_\_\_\_ 3 temperature settings
- \_\_\_\_\_ Illuminated display
- \_\_\_\_\_ Detachable cord
- \_\_\_\_\_ BSS (Beurer Safety System)
- \_\_\_\_\_ Removeable cotton cover - washable
- \_\_\_\_\_ Oeko-Tex 100
- \_\_\_\_\_ The textiles used for this device meet the stringent human ecological requirements of Oeko-Tex Standard 100, as verified by Hohenstein Research Institute.
- \_\_\_\_\_ 100 watt
- \_\_\_\_\_ Product measurements: 40 x 30 cm
- \_\_\_\_\_ Product weight: approx. 460 g

3 year guarantee

Sales unit: 5 / Shipping carton:  
 EAN no.: 4211125 26001 9  
 Item no.: 260.01

