



50 mmol/L


beurer
medical

GL 50 mmol/L Blood glucose monitor

 **german
engineering**



With USB cable and Health Manager software for free download



3 in 1:
Measuring device - Lancing device - Plug in USB



Simple overview of results and trends in target range



Easy-to-use: illuminated display, wide test strips, simple measurement marking enables optimised diabetes management:
before eating, after eating, general marking

Blood volume check: measurement only possible with sufficient sample of blood

Modern design and practical format

White illuminated, easy-to-read LCD display

Lancing device integrated into the unit with adjustable penetration depth

Coloured target range scale & indication on the display to help interpret the measurement result

GlucMemory software in German and English (pre-installed):

Simple overview of results and trends in target range

Health-Manager software download for free under:

www.beurer.com > SERVICE > Download

Starter set: 3-in-1 monitor, AST cap, 10 lancets, 10 individually packed test strips, batteries, USB cover, pouch, instructions for use, brief instructions, blood glucose diary

Suitable for self-testing

Only 0.6 µl of blood necessary

Period of measurement: approx. 5 seconds

Codefree

480 stored results with date, time (data memory also with battery change) and measurement markings

Average value calculation for 7, 14, 30 or 90 days for
- all measurements - marked values

Blood sample: capillary whole blood, venous whole blood

Plasma calibrated

Measurement range: 1.1 - 35.0 mmol/L

Complies with EN ISO 15197:2015 in its entirety (obligatory since Mai 2016)

Automatic switch off

Product measurements (L/H/D): 123 x 28 x 16 mm

Product weight: approx. 36 g (with batteries)



50 mmol/L

beurer
medical

2 x Lithium batteries CR 2032 for more than 1.000 measurements

Storage/Transport:	Temperature:	Rel. humidity:
Unit	-10°C to +65°C	< 90 %
Test stripes	+ 2°C to +30°C	< 90 %
Operating	+ 10°C to +40°C	< 90 %

Medical product, 5 year guarantee

Sales unit: 5 / Shipping carton: 4 x 5

	EAN no.:	Item no.:
Deep Black	4211125 463212	463.21
Pure White	4211125 463250	463.25



Pure white





Deep black



optional



CONNECTION



SOFTWARE

Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via *Bluetooth*®.



Syns with
 • iOS ≥ 9.0
 • Android™ ≥ 5.0
 • Bluetooth® ≥ 4.0



List of supported devices

GL50_mmol_L_2021-02-25_01_DS_EN_BEU Subject to errors and changes

The *Bluetooth*® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.